

# Robotic protocol for RNA extraction by MAVRCIS

## Reagents

1. Binding buffer:
  - 10 mM Bis-Tris /HCl pH< 6.5, 3M guanidium hydrochloride, in 90% Ethanol
2. SiMNP in 96-well plate
3. Wash buffer:
  - Binding buffer + TRIzol at 1:1 v/v
4. Ethanol wash
  - 90% Ethanol

## Materials

1. Magnet for bead separation in plates
  - V&P Scientific, Inc. Cat. VP 771MDWZM-1-ALT
2. Patient Sample plate (96 deepwell)
  - Any deepwell plate
  - Same plates can be used for 2 waste positions
3. Extraction-plate (96 deepwell, square-well, round-bottom)
  - Whatman, Cat. 7701-5200
4. Microplate for bead input (regular 96-well, round-bottom)
  - e.g. Corning Cat. 3788
5. Universal microplate lids
  - Whatman / Cytivia Cat. 77041001
  - Fit all the plates used
6. Evo MCA96 nested tips 200µl sterile
  - Axygen, Cat. EV-200-NTR-S

## Notes

1. Trizol contains phenol and guanidinium isothiocyanate and needs to be discarded accordingly.

## TECAN EVO Script Summary

- volumes are per well
- all plates (except 2<sup>nd</sup> waste position) are covered with lids, to be removed before each pipetting step and placed back afterwards

<i>RNA Binding</i>	Pickup MCA tips
	Add 2 x 150µl binding buffer to extraction plate
	Briefly shake microplate containing beads (resuspend)
	Transfer 40µl beads solution to extraction plate (pipette up / down on aspiration)
	Transfer 2 x 150µl patient sample from input plate to extraction plate
	Shake 5 min @ 1300 rpm
	Move extraction plate to magnet, 30s – 1 min settling time
	Aspirate and discard supernatant
	Move extraction plate back to regular pipetting position
	Change MCA tips
<i>Wash I</i>	Add 2 x 150µl wash buffer
	Shake 20s @ 1300 rpm
	Move extraction plate to magnet, 30s – 1 min settling time
	Aspirate and discard supernatant
	Move extraction plate back to regular pipetting position
	Change MCA tips
<i>Ethanol Wash</i> <b>repeat 4x</b>	Add 2 x 150µl 90% EtOH
	Shake 20s @ 1300 rpm
	Move extraction plate to magnet, 30s – 1 min settling time
	Aspirate and discard supernatant
	Move extraction plate back to regular pipetting position
	Change MCA tips (optional between steps)
<i>Drying</i>	Move extraction plate to heat block, incubate 20 min @ 50 C
	Move extraction plate back to regular pipetting position
<i>Elution</i>	Add 40µl nuclease-free water to extraction plate wells
	Shake 5 min @ 1300 rpm
	Move extraction plate to magnet, 30s – 1 min settling time
	Move destination plate from hotel to pipetting position, remove lid
	Transfer supernatant from extraction plate to destination plate, place back lid
	Discard tips